Sapphire Housing

Note: Document originally drafted in the English language.

The Sapphire housing is designed for Varian Rad series rotating anode inserts having 100mm (4 inch) diameter targets.

Product Description

IEC Classification Class 1, Type B

Weight, Approximate:
Housing & Tube 18.9 kg (41.6 lbs)

Mounting Port plate or Trunnion
Standard Housing - 1/4 X 20 Screws
Metric Housing - M6-1 Screws

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Classification CEI Classe 1, Type B

Poids, Approximatif:
Gaine et Tube 18,9 kg (41.6 lbs)

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Product Description	
Maximum Potential Difference	150 kV
Cathode to Ground	82.5 kV
Anode to Ground	82.5 kV
Maximum X-Ray Tube Assembly Heat Content	1,111 kJ (1,500 kWh)
Maximum Continuous Heat Dissipation	
without air circulator	216 W (300 HU/sec)
with air circulator	432 W (600 HU/sec)
X-Ray Tube Assembly (Insert and Housing):	
Permanent Filtration	0.7mm Al/75kV IEC 60522
Loading Factors for Leakage Radiation	150 kV, 4.0 mA
Thermal Switch	Normally Closed contact rating - 7A @ 120Vac or 30 Vdc Max.
Open	85°C ±3.9°C (185°F ±7°F)
Closed	73.9°C ±2.8°C (165°F ±5°F)
Federal Standard High Voltage Cable	Complies to IEC 60526
Focal Spot to Port Distance	2.2 inches (56 mm)
X-Ray Source Assembly	
Sapphire	Complies to IEC 60601-2-28

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



Filter Assembly
Montage de Filtre
Filteraufbau
Filtro Asamblado

Optional
Filter Kit - IK163
P/N 7576

Note: Window aperture is .875" (22mm) square on 16° mode and .656" (16.7mm) square on 10° and 12° models.

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Stator Power:
Immediately following high speed anode rotation, the rotor speed must be reduced to 3000 r/min or less within 10 seconds using a suitable dynamic braking device.

No more than two high speed starts per minute are permissible. The starting voltage must never exceed 600 volts rms or 1200 volts peak.

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Stator Cord: A six wire shield-grounded stator cord should be connected from the stator terminal to the motor control. All listed voltage and amperage measurements should be taken at the housing end.

If the stator cord is not provided, the cord used must be UL listed or a UL recognized component for all products used in the USA. Outside the USA, the cord used must comply with all applicable regulatory and statutory requirements for electrical and safety. The cord must have the following characteristics: six conductor, 18AWG (0.9mm²) Cu, rated 600V, 90°C minimum, shielded cable with PVC jacket 0.032 inch (0.081mm) thick minimum. Lead terminations are provided and must be used with wire that matches the above description.

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