SMARTLABTOOLS™ QUALITY CONTROL SYSTEM

Parallel Testing and Implementing a New Lot of QC Materials

Templates: SLT_413, SLT_417 SLT_400

Daniel W. Leighton, MS,MT(ASCP),HCLD(ABB),CLB dan@smartlabtools.com

Ver.113016

Objective: Parallel Testing New QC Lot(s)

- SLT Templates are used to facilitate the calculations necessary to establish QC limits for a new lot of controls.
- The examples shown are for QC of the ACCESS-II immunochemistry analyzer using BioRad Controls.
- Note: CLIA regulations require that the laboratory establish it's own mean and standard deviation through repetitive testing. 493.1218 (5d)

Step-1 Run New Set of Controls x 10 (as unknowns)

For cross-over between two different lots of control materials, lab may calculate the mean for the new material from the first 10 measurements, and use the CV from the previous lot of QC material to calculate the SD, then used to calculate QC limits.

New control materials should be run in parallel as unknowns over multiple days/runs

403/1	BIORADLOWNEWLC	Serum	annessa pi posta " e propio de	n di manda dat tentera jan data.
	40891	1	hFSH hLH Ferritin PRL TotT4 FOLW TotT3 PSA-Hyb TSH VitB12 Testo FT3 FRT4 E2	9.38 mIU/mL 4.39 mIU/mL 22.0 ng/mL 7.83 ng/mL 8.59 ug/dL 2.26 ng/mL 1.01 ng/mL 0.32 ng/mL 0.61 uIU/mL 192 pg/mL 0.97 ng/mL 2.15 pg/mL 0.78 ng/dL 51 pg/mL
403/2	BIORADHIGHNEWL(40893	Serum 1	hFSH hLH Ferritin PRL TotT4 FOLW TotT3 PSA-Hyb TSH VitB12 Testo FT3 FRT4 E2	33.48 mIU/mL 65.08 mIU/mL 266.9 ng/mL 41.34 ng/mL 19.13 ug/dL 12.06 ng/mL 2.72 ng/mL 24.96 ng/mL 22.36 uIU/mL 623 pg/mL 10.46 ng/mL 7.90 pg/mL 4.48 ng/dL 747 pg/mL

1 of 10 Measurement Replicates

Step 2. Use SLT 413 Templates to Analyze Data (Level-1 QC Shown)

Smart Smart	PR	ECISION	AND AC	CURACY	STATIS	TICAL AS	SESSME	NT				
Cabiconstructor			INTERNA		INE ASS	OCIATES						
?	ACCES	S-II PARA	LLEL TEST	TING NEW	LOT OF B	IORAD IA-I	PLUS CON	TROLS				
Analyte :	FSH	LH	FERRITIN	PROLACT	TT4	FOLATE	Π3	PSA				
QC Material :	BIORAD											
Lot Number :	40891											
Expiration :	6/17											
Target Values :	8.2	4.27	19.7	6.94	7.99	2.75	0.862	0.33				
Ranges :	6.45-9.96	3.36-5.19	15.7-23.7	5.80-8.08	6.09-9.89	1.18-4.32	0.398-1.33	1.82-2.74				
Run	L-1	L-1	L-1	L-1	L-1	L-1	L-1	L-1				
1	8.78	4.31	24.1	7.67	8.86	2.23	1.00	0.33				
2	9.38	4.39	22.0	7.83	8.59	2.26	1.01	0.32				
3	9.13	4.29	22.4	7.90	8.47	2.27	1.00	0.31				
4	9.50	4.47	23.6	7.85	8.59	2.02	0.97	0.33				
5	8.17	4.03	20.9	7.79	8.99	2.50	0.95	0.33				
6	8.52	4.58	20.5	8.40	8.58	2.46	0.96	0.34				
7	7.70	4.45	21.0	8.13	8.90	2.60	0.95	0.36				
8	7.30	4.18	23.8	8.36	9.07	2.67	1.01	0.36				
9	8.40	3.78	19.2	7.85	9.78	2.59	0.92	0.34				
10	7.25	3.04	19.3	7.43	8.90	2.39	0.86	0.31				
11												
12												
13	N/0	0.000	f		r - t		d fo					
14		ans	пог	n m	sen	Use	a 10					
15	(T -		" \/_	l								
16	61	irget	va	lues								
17												
18												
19												
20												
N :	10	10	10	10	10	10	10	10				
Mean :	8.41	4.15	21.68	7.92	8.87	2.40	0.96	0.33				
1 SD :	0.81	0.46	1.79	0.30	0.38	0.20	0.05	0.02				
% CV :	9.65	10.96	8.27	3.78	4.26	8.48	4.90	5.31				
Target Value :	8.20	4.27	19.70	6.94	7.99	2.75	0.86	0.33				
% Recovery :	102.60	97.24	110.05	114.14	111.05	87.24	111.72	100.91				
Clear Form	Reset	Reset	Reset	Reset	Reset	Reset	Reset	Reset				
Comments : PAR	Comments : PARALLEL FOR NEW LOT OF CONTROLS VS. E-INSERT											
Analyst : BRITTA	NY		t	1/22/2016	Approved	by : DL						

Smart	PR	PRECISION AND ACCURACY STATISTICAL ASSESSMENT													
CarplectorremateSingle			INTERNA			OCIATES									
?	ACCES	S-II PARA	LLEL TEST	TING NEW	LOT OF B	IORAD IA-	PLUS CON	TROLS							
Analyte :	TSH	VIT-B12	TESTO	FT3	FT4	E2									
QC Material :	BIORAD														
Lot Number :	40891														
Expiration :	6/17														
Target Values :	0.732	196	0.875	2.28	0.773	40.3									
Ranges :	0.539-0.925	143-249	0.647-1.10	1.82-2.74	0.598-0.949	<20-84.8									
Run	L4	L-1	L-1	L-1	L-1	L-1	L-1	L-1							
1	0.61	192	0.97	2.15	0.78	51									
2	0.60	195	0.93	2.19	0.80	50									
3	0.66	183	0.94	2.25	0.82	44									
4	0.70	166	0.94	1.94	0.82	48									
5	0.64	187	0.91	2.05	0.78	61									
6	0.66	190	0.92	2.22	0.92	53									
7	0.61	190	1.01	2.29	0.85	48									
8	0.60	203	1.06	2.22	0.89	57									
9	0.61	191	1.07	2.19	0.85	42									
10	0.64	182	0.89	2.27	0.83	54									
11															
12															
13	N I - 4	FC													
14	NOt	e Eż	. %h	(eco	verv	/ IS									
15															
16	Hia	1. bt	it wi	thin	"Ins	ert"	Limi	ts							
17	9	.,													
18															
19															
20															
N :	10	10	10	10	10	10									
Mean :	0.63	187.90	0.96	2.18	0.83	50.80									
1 SD :	0.03	9.73	0.06	0.11	0.05	5.75									
% CV :	5.22	5.18	6.49	4.93	5.43	11.32									
Target Value :	0.73	196.00	0.88	2.28	0.77	40.30	/								
% Recovery :	86.48	95.87	110.17	95.48	107.89	126.05									
Clear Form	Reset	Reset	Reset	Reset	Reset	Reset	Reset	Reset							
Comments : PAF	RALLEL TESTI	NG NEW LOT	OF CONTROL	S VS. E-INSE	RT										
Analyst : BRITT	ANY		1	1/22/2016	Approved	Approved by : D. LEIGHTON									

Step 2. (cont.) Use SLT 413 Templates to Analyze Data (Level-3 QC Shown)

Smart Smart	PR	ECISION	AND AC	CURACY	STATIS	TICAL AS	SESSME	NT
Complex forms made Simple			INTERNA	L MEDIC	INE ASS	OCIATES		
	ACCES	S-II PARA	LLEL TEST	TING NEW	LOT OF B	IORAD IA-	PLUS CON	ITROLS
Analyte :	FSH	LH	FERRITIN	PROLACT	TT4	FOLATE	Π3	PSA
QC Material :	BIORAD							
Lot Number :	40893							
Expiration :	6/17							
Target Values :	31.1	55.3	269	39.0	18.6	14.9	2.48	25.1
Ranges :	24.6-37.6	44.6-65.9	174-364	33.0-44.9	15.3-22.0	7.55-22.2	1.59-3.36	19.1-31.1
Run	L-3	L-3	L-3	L-3	L-3	L-3	L-3	L-3
1	33.48	65.08	266.9	41.34	19.13	12.06	2.72	24.96
2	34.74	61.51	271.1	41.66	19.21	10.85	2.64	25.92
3	35.48	64.45	292.9	42.28	19.73	11.85	2.73	24.23
4	34.44	63.74	306.8	41.07	18.15	11.64	2.70	25.04
5	30.90	57.87	305.9	40.91	19.59	12.27	2.68	24.97
6	33.77	60.40	308.4	44.24	18.00	13.22	2.64	26.10
7	33.88	60.29	288.2	43.99	19.08	13.19	2.67	26.46
8	33.23	51.33	338.6	43.15	18.94	12.91	2.72	25.60
9	33.44 54.53		282.8	42.70	19.08	12.89	2.78	15.37
10	28.46	48.85	277.7	41.54	19.79	12.83	2.60	23.86
11								
12								
13								
14								
15								
16								
17								
18								
10								
20								
20	10	10	10	10	10	10	10	10
Moan :	10	50.01	202.02	43.20	10.07	10 12 27	2.60	24.25
Medil :	33.10	50.01	233.93	42.29	19.07	0.77	2.09	24.20
TSD:	1 SD: 2.05		21.57	1.20	0.60	0.77	0.05	3.22
% CV :	6.18	9.52	1.34	2.83	3.15	6.25	1.97	13.30
Target Value :	31.10	55.30	269.00	39.00	18.60	14.90	2.48	25.10
% Recovery :	106.69	106.34	109.27	108.43	102.53	83.03	108.39	96.62
Clear Form	Reset	Reset	Reset	Reset	Reset	Reset	Reset	Reset
Comments : PAR	ALLEL TEST	NG NEW LOT	OF CONTROL	S VS. E-INSE	RT			
Analyst : BRITTA	NY		1	1/22/2016	Approved	by : D. LEIGH	TON	

Smart	PRECISION AND ACCURACY STATISTICAL ASSESSMENT												
Complex forms made Simple			INTERN/	AL MEDIC	INE ASS	OCIATES							
?	ACCES	S-II PARA	LLEL TES	TING NEW	LOT OF B	IORAD IA-I	PLUS CON	TROLS					
Analyte :	TSH	VIT B-12	TESTO	FT3	FT4	E2							
QC Material :	BIORAD												
Lot Number :	40891												
Expiration :	6/17												
Target Values :	26.1 573		9.70	7.71	4.29	798							
Ranges :	19.7-32.4	412-733	7.66-11.7	6.17-9.25	3.45-5.14	544-1051							
Run	L-3	L-3	L-3	L-3	L-3	L-3	L-3	L-3					
1	22.36	623	10.46	7.90	4.48	747							
2	25.00	611	9.91	7.65	4.25	764							
3	19.36	613	10.11	7.83	4.26	813							
4	22.01	668	10.32	8.23	4.35	806							
5	21.10	599	9.73	7.53	4.14	798							
6	20.82	615	9.96	8.05	4.98	824							
7	22.30	652	10.52	7.92	4.64	858							
8	23.78	609	10.36	8.09	4.47	748							
9	20.61	634	10.23	7.75	4.43	805							
10	23.97 575		10.29	8.02	4.42	733							
11													
12													
13													
14													
15													
16													
17													
18													
19													
20													
N :	10	10	10	10	10	10							
Mean :	22.13	619.90	10.19	7.90	4.44	789.60							
1 SD :	1.74	26.40	0.26	0.21	0.24	39.99							
% CV :	7.86	4.26	2.51	2.70	5.31	5.06							
Target Value :	26.10 573.00		9.70	7.71	4.29	798.00							
% Recovery :	84.79 108.19		105.04	102.43	103.54	98.95							
Clear Form	Reset	Reset	Reset	Reset	Reset Reset Reset Reset								
Comments : PAR	Comments : PARALLEL TESTING NEW LOT OF CONTROLS VS. E-INSERT												
Analyst : BRITTA	NY		1	11/22/2016	Approved by : D. LEIGHTON								

Step 3. Compare Preliminary Lab Values vs. Insert .. *Mean & Range are used on the SLT_413 Form*

Means from parallel study should fall within the manufacturer's stated range.

Insert limits should be used only as guides in setting initial control limits for testing new control materials.

Liquichek[™] Immunoassay Plus Control Levels 1, 2 and 3



http://www.myeinserts.com/40890

Revision Date 2016-11-17 → Indicates Revised Information

INSTRUMENT (1)

BIO-RAD

		L	evel 1 - 40891	Le	evel 2 - 40892	Le	vel 3 - 40893
	Units	Mean	Range	Mean	Range	Mean	Range
BECKMAN COULTER ACCESS / 2 / 2i							
Estradiol	pg/mL	40.3	<20.0 - 84.8	339	231 - 447	798	544 - 1051
Ferritin	ng/mL	19.7	15.7 - 23.7	121	93.4 - 148	269	174 - 364
Folate (FOLW)	ng/mL	2.75	1.18 - 4.32	10.2	5.50 - 14.9	14.9	7.55 - 22.2
Follicle Stimulating Hormone (FSH) (hFSH)	mlU/mL	8.20	6.45 - 9.96	18.7	14.7 - 22.6	31.1	24.6 - 37.6
Luteinizing Hormone (LH) (hLH)	mlU/mL	4.27	3.36 - 5.19	18.7	14.3 - 23.1	55.3	44.6 - 65.9
Prolactin	ng/mL	6.94	5.80 - 8.08	16.0	13.3 - 18.6	39.0	33.0 - 44.9
PSA (Total) (Hybritech PSA)	ng/mL	0.330	0.250 - 0.410	3.71	2.82 - 4.60	25.1	19.1 - 31.1
T3 (Free) (Free T3)	pg/mL	2.28	1.82 - 2.74	5.37	4.25 - 6.45	7.71	6.17 - 9.25
T3 (Total) (Total T3)	ng/mL	0.862	0.398 - 1.33	1.72	1.05 - 2.38	2.48	1.59 - 3.36
T4 (Free) (Free T4)	ng/dL	0.773	0.598 - 0.949	2.56	2.07 - 3.06	4.29	3.45 - 5.14
T4 (Total) (Total T4)	µg/dL	7.99	6.09 - 9.89	12.2	9.61 - 14.7	18.6	15.3 - 22.0
Testosterone	ng/mL	0.875	0.647 - 1.10	4.57	3.52 - 5.62	9.70	7.66 - 11.7
Thyroid Stimulating Hormone (TSH) (hTSH, HYPERsensitive)	µIU/mL	0.732	0.539 - 0.925	5.07	4.05 - 6.09	26.1	19.7 – 32.4

6

Step 4. Compare with Manufacturer's Peer Report

Unity Manufacturer Report for Beckman Coulter

Immunoassay Plus • Lot 40890 • Exp 30-Jun-2017

Estradiol, E	2 Chem	ilumineso	ence po	g/mL						
	Level	Mon	Cum	Level	l Mon	Cum	Level	Mon	Cum	
Beckman Co	oulter Ac	cess, LXi	725, Dx	C 600i IA System	s		_			
Mear	14	51.06	52.08	9	376.6	373.3	2	825.3	818.0	
SE)	12.76	12.75	4	38.62	33.04	ں ا	59.04	60.89	
C/	/	25.0	24.5		10.3	8.8		7.2	7.4	
# Points	5	389	3332		152	1454		355	3195	
# Labs	5	15	21		8	11		14	18	

The lab will use the calculated mean, along with the peer CV to calculate it's interim QC limits for the new QC materials.

We note that the peer Level-1 E2 value more closely matches the E2 value from the lab study (50.8). The most current peer reports contain the more reliable comparative values.

Step 5. Use Template SLT 417 to Calculate 2SD Limits

Lab determined Means, and Peer CV's are used here to calculate new 2SD Limits.

Optionally, Lab Previous CV's are used.(HCV)

QC Limits are best determined by cumulative statistics from 3-6 months testing.

Smart LabTools		(CALC	ULA	TE Q	CLIM	ITS L	JSIN	g his	TORIC	CALC	:V% (HCV)	
Complex Forms made Simple	?			IN	TERN	AL M	EDIC	INE	ASSO	CIAT	ES LA	В			
METHOD :	ACCES	S-2									ACCES	S-2			
CONTROLS :	BIORAD		-1								BIORAD LEVEL-3				
LOT #'s :	40891										40893				
EXPIRATION :	6/17		_								6/17				
ANALYTE	MEAN HCV SD -2SD +2SD MEAN HCV SD -2SD						-2SD	+2SD	MEAN	HCV	SD	-2SD	+2SD		
ESTRADIOL-2	50.8	24.5	12.45	25.91	75.69						790	7.4	58.46	673.08	906.92
FERRITIN	21.7	6.9	1.50	18.71	24.69						294	6.9	20.29	253.43	334.57
FOLATE	2.40	8.3	0.20	2.00	2.80						12.4	7.1	0.88	10.64	14.16
FSH	8.41	6.5	0.55	7.32	9.50						33.2	8.8	2.92	27.36	39.04
FT4	0.83	6.1	0.05	0.73	0.93						4.29	4.4	0.19	3.91	4.67
LH	4.15	7.7	0.32	3.51	4.79						58.8	5.4	3.18	52.45	65.15
PROLACTIN	7.92	5.2	0.41	7.10	8.74						42.3	4.5	1.90	38.49	46.11
PSA, HYB	0.33	6.6	0.022	0.286	0.374						24.25	4.6	1.116	22.019	26.481
тѕн	0.63	5.6	0.04	0.56	0.70						22.13	6.1	1.35	19.43	24.83
ттз	0.96	9.6	0.09	0.78	1.14						2.69	5.9	0.16	2.37	3.01
ТТ4	8.87	7.8	0.69	7.49	10.25						19.1	6.1	1.17	16.77	21.43
VIT B12	188	8.8	16.54	154.91	221.09						620	6.9	42.78	534.44	705.56
FREE T3	2.18	7.3	0.16	1.86	2.50						7.90	5.8	0.46	6.98	8.82
TESTOST	0.96	8.7	0.08	0.79	1.13						10.2	5.9	0.60	9.00	11.40
Means from															
initial Parallel															
study 11/2016															
CV is peer															
of 10/2016															

Step 6. Enter 2SD Limits into QC Assessment Template

When using the SLT Daily QC Assessment Template (SLT_400), 2SD Limits are Entered on the 2nd page of the 2-page Template.

Smart	NTE	RNA	L ME	DICI	NE A	ASSC		TES	LAB	7	QC L	esson
Complex Forms made Simple ES	TABL	ISHE	D QU	ALIT	Y CO	NTRC	DL PA	RAM	ETER	s 🔍	Rese	et All
TEST SYSTEM:	BECKN	IAN AC	CESS 2	2	BECKMAN ACCESS 2				PEER MEAN/HCV			
CONTROLS:	BIORA	D IA-PL	US LEV	/EL-1	BIORAD IA-PLUS LEVEL-3				AS TARGET/QC LIMITS			
LOT NUMBERS:	40891				40893							
EXPIRATION:	6/30/17	7			6/30/1	7						
Analyte Description	L-1 -2SD	L-1 +2SD	L-1 Mean	L-1 1SD	L-2 -2SD	L-2 +2SD	L-2 Mean	L-2 1SD	L-3 -2SD	L-3 +2SD	L-3 Mean	L-3 1SD
ESTRADIOL 2	25.9	75.7	50.80	12.45	673	907	790.00	58.50				
FERRITIN	18.7	24.7	21.70	1.50	253	335	294.00	20.50				
FOLATE	2.0	2.8	2.40	0.20	10.6	14.2	12.40	0.90				
FSH	7.32	9.50	8.41	0.55	27.3	39.0	33.15	2.93				
FT4	0.73	0.93	0.83	0.05	3.91	4.67	4.29	0.19				
LH	3.51	4.79	4.15	0.32	52.4	65.2	58.80	3.20				
PROLACTIN	7.1	8.74	7.92	0.41	38.5	46.1	42.30	1.90				
PSA, HYB	0.286	0.374	0.33	0.02	22.0	26.5	24.25	1.13				
TSH	0.56	0.70	0.63	0.04	19.4	24.8	22.10	1.35				
ТТ3	0.78	1.14	0.96	0.09	2.37	3.01	2.69	0.16				
TT4	7.49	10.25	8.87	0.69	16.77	21.43	19.10	1.17				
VIT B12	155	221	188.00	16.50	534	706	620.00	43.00				
FREE T3	1.86	2.50	2.18	0.16	6.98	8.82	7.90	0.46				
TESTOST	0.79	1.13	0.96	0.09	9.0	11.40	10.20	0.60				

Step 7. Daily QC Statistical Assessment (SLT_400)

QC Form (page-1) is now ready to test data from Study to rule out typo's.

Verified QC Limits are set into Analyzer & LIS QC programs.

	?	I	NTE D/		AL ′Q.	MEI c. si	MEDICINE ASSOCIATES LAB									?			
TEST SYSTEM:	BECK	MAN AC	CESS	2		BECK	MAN AG	CESS	2		PEER MEAN/HCV					Bias # CTLs			
CONTROLS:	BIORA	D IA-PL	US LE	VEL	1	BIORAD IA-PLUS LEVEL-3				AS TARGET/QC LIMITS				2					
LOT NUMBERS:	40891					40893								Trend Fla	ag =				
EXPIRATION:	6/30/1	7				6/30/17									1.5				
Analyte Description	L-1 Mean	Test Value	Bias	SDI (Z)	QC In?	L-2 Mean	Test Value	Bias	SDI (Z)	QC In?	L-3 Mean	Test Value	Bias	SDI (Z)	QC In?	Ave SDI (Z)	Trend Alert		
ESTRADIOL 2	50.80	50	-0.80	-0.06	In	790.00	764	-26.00	-0.44	In						-0.25			
FERRITIN	21.70	24.1	2.40	1.60	In	294.00	271	-23.00	-1.12	In						0.24	*		
FOLATE	2.40	2.23	-0.17	-0.85	In	12.40	10.9	-1.50	-1.67	In						-1.26	*		
FSH	8.41	8.78	0.37	0.68	In	33.15	34.7	1.55	0.53	In						0.60			
FT4	0.83	0.80	-0.03	-0.60	In	4.29	4.25	-0.04	-0.21	In						-0.41			
LH	4.15	4.31	0.16	0.50	In	58.80	61.5	2.70	0.84	In						0.67			
PROLACTIN	7.92	7.67	-0.25	-0.61	In	42.30	41.7	-0.60	-0.32	In						-0.46			
PSA, HYB	0.33	0.33	0.00	0.00	In	24.25	25.9	1.65	1.47	In						0.73			
тѕн	0.63	0.60	-0.03	-0.86	In	22.10	22.4	0.30	0.22	In						-0.32			
ТТ3	0.96	1.00	0.04	0.44	In	2.69	2.64	-0.05	-0.31	In						0.07			
TT4	8.87	8.86	-0.01	-0.01	In	19.10	19.2	0.10	0.09	In						0.04			
VIT B12	188.00	195	7.00	0.42	In	620.00	611	-9.00	-0.21	In						0.11			
FREE T3	2.18	2.19	0.01	0.06	In	7.90	7.65	-0.25	-0.54	In						-0.24			
TESTOST	0.96	0.93	-0.03	-0.35	In	10.20	9.91	-0.29	-0.48	In						-0.42			

Step 8. Inform Analysts Of Lot Change

- Post Clear Messages as to which QC Materials, QC Files, QC Templates are to be used once change is made.
- AVOID mix-ups such as using old material with new QC Files / Templates.. or... new QC material with old QC Files / Templates.
- BEST REMOVE OLD QC MATERIALS to avoid such incidents & resultant headaches.
- Closely monitor QC activities for first few days new QC materials are in use.
- Make on-going statistical parameter adjustments if warranted.

